

Report No. 24PTP_01PM04240324

Page 1 of 3

Preventive Maintenance Action Report

EQUIPMENT : VASCULAR SCREENING UNIT (ABI)

ID CODE : PLP03861

MANUFACTURER : FUKUDA DENSHI

MODEL : VASERA VS-2000

SERIAL No. : 50300417

LOCATION : General Surgery

SUBMITTED BY : PHYATHAI PHAHOLYOTHIN HOSPITAL

670/1 Phaholyothin Road, Samsen Nai, Phaya Thai, Bangkok 10400

Tel: 02-271-7000 Fax: -

PM. DATE : 25 APRIL 2024

ISSUE DATE : 1 MAY 2024

Performed by: Approved by:

USANEE CHALONGPHOKSILCHAI PHUWASIN YIWKIM

CONDITION OF THIS RESULT OF TEST

1. REFERENCE STANDARD INSTRUMENT:

MASTERMANUFACTURERMODELSERIAL NO.CERTIFICATE NO.DUE DATEElectrical SafetyFlukeESA615349403423MD7747 Jun 2024

Analyzer

2. THIS CERTIFICATION IS TRACEABLE TO:

- Technology Promotion Association (Thai-Japan)

3. THIS RESULT OF TEST WAS FOUND ACCURATE AS SHOW ON DATE AND PLACE OF TEST ONLY



Report No. 24PTP_01PM04240324

Page 2 of 3

Preventive Maintenance Action Report

EQUIPMENT : VASCULAR SCREENING UNIT (ABI)

ID CODE : PLP03861 SERIAL No. : 50300417

1. Qualitative Tests

No.	List	Pass	Fail	None	Comment
1	Chassis/Housing	/			
2	Mount/Fasteners	/			
3	Casters/Brakes	/			
4	AC Plug	/			
5	Line Cord	/			
6	Strain Reliefs	/			
7	Circuit Breaker/Fuse	/			
8	Tubes/Hoses	/			
9	Fittings/Connectors	/			
10	Transducers	/			
11	Controls/Switches	/			
12	Battery/Charger	/			
13	Time/Date Settings	/			
14	User Calibration/Self-Test	/			
15	Alarms	/			
16	Labeling	/			
17	Accessories	/			

2. Quantitative Tests

No.	List	Pass	Fail	None	Comment
1	Grounding Resistance ≤0.5 Ω	/			0.103 Ω
2	Chassis Leakage Current ≤500 μA	/			64.7 μΑ
3	Lead-to-Ground Leakage Current ≤100 μA			/	
	grounded; ≤500 μA ungrounded				
4	Air Leakage ≤15 mmHg/min	/			
5	Heart Rate ±10%	/			
6	Static Pressure Accuracy ±3 mmHg or	/			
	*Oscillometric Pressure ±8 mmHg				
7	Maximum Pressure ≤330 mmHg	/			

3. Preventive Maintenance

No.	List	Pass	Fail	None	Comment
1	Clean the exterior and accessories	/			



Report No. 24PTP_01PM04240324

Page 3 of 3

Preventive Maintenance Action Report

EQUIPMENT : VASCULAR SCREENING UNIT (ABI)

ID CODE : PLP03861 SERIAL No. : 50300417

3. Preventive Maintenance (Continues)

No.	List	Pass	Fail	None	Comment
2	Lubricate			/	
3	Adjust			/	
4	Replace			/	

Electrical Safety Test: Master used Electrical Safety Analyzer (Metron / QA-90) according to IEC.60601.1