

Report No.

24PYT3_DECPM12240110

PUKARIN TONGKLIANG

Page 1 of 2

Preventive Maintenance Action Report

EQUIPMENT	: ASPIRATORS, EMERGENCY
ID CODE	: PYT3_00138
MANUFACTURER	: MEDELA
MODEL	: VARIO 18 C/I
SERIAL No.	: 1212016
LOCATION	: OPERATING ROOM
SUBMITTED BY	: PHYATHAI 3 HOSPITAL
	111 Phet Kasem Rd., Pak Khlong Phasi Charoen, Phasi Charoen Bangkok 1016
	Tel: (662) 467-1111 Fax: (622) 467-1111
PM. DATE	: 3 DECEMBER 2024
ISSUE DATE	: 30 DECEMBER 2024
Performed by:	Approved by :
i cii ciiiica c y .	1 ipprovou by .

CONDITION OF THIS RESULT OF TEST

1. REFERENCE STANDARD INSTRUMENT:

SOM KONKAEW

MASTER MANUFACTURER MODEL SERIAL NO. CERTIFICATE NO. DUE DATE

2. THIS CERTIFICATION IS TRACEABLE TO:

3. THIS RESULT OF TEST WAS FOUND ACCURATE AS SHOW ON DATE AND PLACE OF TEST ONLY



Report No. 24PYT3_DECPM12240110

Page 2 of 2

Preventive Maintenance Action Report

EQUIPMENT : ASPIRATORS, EMERGENCY

ID CODE : PYT3_00138

SERIAL No. : 1212016

1. Qualitative Tests

No.	List	Pass	Fail	None	Comment
1	Chassis/Housing	/			
2	Tubes/Hoses	/			
3	Fittings/Connectors	/			
4	Filters			/	
5	Controls/Switches	/			
6	Overflow Protection			/	
7	Intermittent Operation			/	
8	Indicators/Display	/			
9	Labeling	/			
10	Accessories	/			

2. Quantitative Tests

No.	List	Pass	Fail	None	Comment
1	Maximum Free Flow surgical/tracheal ≥25 L/min,			/	
	thoracic ≥20 L/min				
2	Maximum Vacuum Low-volume>100 mmHg;			/	
	thoracic >50 mmHg ; surgical/tracheal ≥400 mmHg				
3	Vacuum Gauge Accuracy ±5% FSO(Full Scale	/			
	output)				

3. Preventive Maintenance

No.	List	Pass	Fail	None	Comment
1	Clean the exterior and accessories	/			
2	Lubricate			/	
3	Adjust			/	
4	Replace			/	

4. Remark Chassis Leakage Current ≤500 μA : 76.7 μA