

Page 1 of 4

Preventive Maintenance Action Report

EQUIPMENT : SCANNING SYSTEMS, ULTRASONIC, GENERAL-PURPOSE

ID CODE : PTP01740

MANUFACTURER : GE HEALTHCARE

MODEL : LOGIQ C5 PREMIUM

SERIAL No. : 290139WX5

LOCATION : Obstetric & Gynecology Center

SUBMITTED BY : PHYATHAI PHAHOLYOTHIN HOSPITAL

670/1 Phaholyothin Road, Samsen Nai, Phaya Thai, Bangkok 10400

Tel: 02-271-7000 Fax: -

PM. DATE : 25 APRIL 2024

ISSUE DATE : 1 MAY 2024

Supatcha

Performed by : _____ Approved by : ____

SUPATCHA BUAPUEAN PHUWASIN YIWKIM

CONDITION OF THIS RESULT OF TEST

1. REFERENCE STANDARD INSTRUMENT:

MASTERMANUFACTURERMODELSERIAL NO.CERTIFICATE NO.DUE DATEMULTI-SUN NUCLEARSONO403802259-5417-1CALPM0923003819 Sep 2024

PURPOSE, MULTI-TISSUE ULTRASOUND PHANTOM

Electrical Safety Fluke ESA615 3494034 23MD774 7 Jun 2024

Analyzer

2. THIS CERTIFICATION IS TRACEABLE TO:

- National Healthcare Systems
- Technology Promotion Association (Thai-Japan)

3. THIS RESULT OF TEST WAS FOUND ACCURATE AS SHOW ON DATE AND PLACE OF TEST ONLY



Page 2 of 4

Preventive Maintenance Action Report

EQUIPMENT : SCANNING SYSTEMS, ULTRASONIC, GENERAL-PURPOSE

ID CODE : PTP01740

SERIAL No. : 290139WX5

1. Qualitative Tests

No.	List	Pass	Fail	None	Comment
1	Chassis/Housing	/			
2	Mount/Fasteners	/			
3	Casters/Brakes	/			
4	AC Plug/Receptacles	/			
5	Line Cord	/			
6	Strain Reliefs	/			
7	Circuit Breaker/Fuse	/			
8	Cables	/			
9	Fittings/Connectors	/			
10	Transducers	/			
11	Filters	/			
12	Controls/Switches	/			
13	Fans	/			
14	Indicators/Displays	/			
15	Time/Date Settings	/			
16	Self-Test	/			
17	Cybersecrity	/			
18	User Calibration	/			
19	Audible Signals	/			
20	Labeling	/			
21	Accessories	/			

2. Mains Voltage

No.	List	Pass	Fail	None	Comment
1	Line (Live) to Neutral ≈ 230 Vac	/			225.7 V
2	Line (Live) to Earth ≈ 230 Vac	/			0.7 V
3	Neutral to Earth ≈ 0 Vac	/			225.6 V

3. Quantitative Tests

No.	List	Pass	Fail	None	Comment
1	Ground Resistance ≤0.5 Ω	/			0.346Ω
2	Chassis Leakage Current ≤ 500 µA [ungrounded]	/			1.7 μΑ
	patient-care equipment				



Page 3 of 4

Preventive Maintenance Action Report

EQUIPMENT : SCANNING SYSTEMS, ULTRASONIC, GENERAL-PURPOSE

ID CODE : PTP01740

SERIAL No. : 290139WX5

3. Quantitative Tests (Continues)

No.	List	Pass	Fail	None	Comment
3	Lead-to-Ground Leakage Current ≤100 μA	/			0.9 μA Probe 1 Name :
	grounded; ≤500 μA ungrounded				C2-5-RC SN:
					297397WX2
4	Lead-to-Ground Leakage Current ≤100 μA	/			0.8 μA Probe 2 Name :
	grounded; ≤500 μA ungrounded				E8CS-RC SN:
					906966WX7
5	Rate Calibration (ECG) 60,80,120 BPM ±5%			/	

4. Probe Test Probe 1

No.	List	Pass	Fail	None	Comment
1	Dead Zone	/			0.38 cm
2	Resolution	/			
3	Measure Caliper Accuracy - Vertical Plane. ± 5%	/			- Vertical Plane. 4 cm.
					= 4.04 cm
4	Measure Caliper Accuracy - Horizontal Plane. ±	/			- Horizontal Plane. 3
	5%				cm. = 3.11 cm
5	Image Uniformity	/			
6	Measure Image Penetration	/			
7	Check Anechoic Image	/			

5. Probe Test Probe 2

No.	List	Pass	Fail	None	Comment
1	Dead Zone	/			0.45 cm
2	Resolution	/			
3	Measure Caliper Accuracy - Vertical Plane. ± 5%	/			- Vertical Plane. 4 cm.
					= 4.01 cm
4	Measure Caliper Accuracy - Horizontal Plane. ±	/			- Horizontal Plane. 3
	5%				cm. = 3.01 cm
5	Image Uniformity	/			
6	Measure Image Penetration	/			_
7	Check Anechoic Image	/			



Page 4 of 4

Preventive Maintenance Action Report

EQUIPMENT : SCANNING SYSTEMS, ULTRASONIC, GENERAL-PURPOSE

ID CODE : PTP01740

SERIAL No. : 290139WX5

8. Preventive Maintenance

No.	List	Pass	Fail	None	Comment
1	Clean the exterior and accessories	/			
2	Lubricate			/	
3	Adjust			/	
4	Replace			/	

Electrical Safety Test: Master used Electrical Safety Analyzer (Metron / QA-90) according to IEC.60601.1

FI-BME-NHS-CP-001/3 Rev.09 Page 4/4 Issued Date 20/07/2023