

Report No. 2

24PYT3\_DECPM12240017

Page 1 of 3

# **Preventive Maintenance Action Report**

| ·              | SOM KONKAEW   | PUKARIN TONGKLIANG   |
|----------------|---------------|--|
| Performed by : | Gen           | Approved by :  |
| ISSUE DATE     | : 30 DECEM    | MBER 2024  |
| PM. DATE       | : 1 DECEME    | BER 2024   |
|                | Tel: (662)    | 467-1111 Fax: (622) 467-1111                                 |
|                | 111 Phet Kase | em Rd., Pak Khlong Phasi Charoen, Phasi Charoen Bangkok 1010 |
| SUBMITTED BY   | : РНҮАТНА     | AI 3 HOSPITAL  |
| LOCATION       | : CHECK-U     | P (ตรวจสุขภาพ)   |
| SERIAL No.     | : JA011162    |  |
| MODEL          | : 53000       |  |
| MANUFACTURER   | : WELCH A     | LLYN   |
| ID CODE        | : PYT3_0139   | 92   |
| EQUIPMENT      | : MONITOR     | RS, PHYSIOLOGIC, VITAL SIGNS                                 |

#### **CONDITION OF THIS RESULT OF TEST**

1. REFERENCE STANDARD INSTRUMENT:

MASTER MANUFACTURER MODEL SERIAL NO. CERTIFICATE NO. DUE DATE

2. THIS CERTIFICATION IS TRACEABLE TO:

3. THIS RESULT OF TEST WAS FOUND ACCURATE AS SHOW ON DATE AND PLACE OF TEST ONLY



*Report No.* 24PYT3\_DECPM12240017

Page 2 of 3

# **Preventive Maintenance Action Report**

EQUIPMENT : MONITORS, PHYSIOLOGIC, VITAL SIGNS

ID CODE : PYT3\_01392

SERIAL No. : JA011162

### 1. Qualitative Tests

| No. | List                       | Pass | Fail | None | Comment |
|-----|----------------------------|------|------|------|---------|
| 1   | Chassis/Housing            | /    |      |      |         |
| 2   | Mount/Fasteners            | /    |      |      |         |
| 3   | Casters/Brakes             | /    |      |      |         |
| 4   | AC Plug                    | /    |      |      |         |
| 5   | Line Cord                  | /    |      |      |         |
| 6   | Strain Reliefs             | /    |      |      |         |
| 7   | Circuit Breaker/Fuse       | /    |      |      |         |
| 8   | Tubes/Bulbs                | /    |      |      |         |
| 9   | Fittings/Connectors        | /    |      |      |         |
| 10  | Transducers                | /    |      |      |         |
| 11  | Controls/Switches          | /    |      |      |         |
| 12  | Battery/Charger            | /    |      |      |         |
| 13  | Time/Date Settings         | /    |      |      |         |
| 14  | Cybersecurity              |      |      | /    |         |
| 15  | User Calibration/Self-Test |      |      | /    |         |
| 16  | Alarms                     | /    |      |      |         |
| 17  | Labeling                   | /    |      |      |         |
| 18  | Accessories                | /    |      |      |         |

## 2. Quantitative Tests

| No. | List                                | Pass | Fail | None | Comment |
|-----|-------------------------------------|------|------|------|---------|
| 1   | Grounding Resistance ≤0.5 Ω         |      |      | /    |         |
| 2   | Chassis Leakage Current ≤500 µA     | /    |      |      | 1.4 μΑ  |
| 3   | Air Leakage ≤15 mmHg/min            | /    |      |      |         |
| 4   | Static Pressure Accuracy ±3 mmHg or | /    |      |      |         |
|     | *Oscillometric Pressure ±8 mmHg     |      |      |      |         |
| 5   | Maximum Pressure (NIBP) ≤330 mmHg   | /    |      |      |         |
| 6   | Heart Rate ±10%                     | /    |      |      |         |

#### 3. Preventive Maintenance

| No. | List                               | Pass | Fail | None | Comment |
|-----|------------------------------------|------|------|------|---------|
| 1   | Clean the exterior and accessories | /    |      |      |         |
| 2   | Lubricate                          |      |      | /    |         |



*Report No.* 24PYT3\_DECPM12240017

Page 3 of 3

# **Preventive Maintenance Action Report**

EQUIPMENT : MONITORS, PHYSIOLOGIC, VITAL SIGNS

ID CODE : PYT3\_01392

SERIAL No. : JA011162

### 3. Preventive Maintenance (Continues)

| No. | List    | Pass | Fail | None | Comment |
|-----|---------|------|------|------|---------|
| 3   | Adjust  |      |      | /    |         |
| 4   | Replace |      |      | /    |         |

FI-BME-NHS-CP-007/3 Rev.04 Page 3/3 Issued Date 20/07/2024